

14" LIQUID NITROGEN FREEZER RACKS

COMPLETE SYSTEMS AND INDIVIDUAL RACKS FOR 14 INCH DIAMETER CONTAINERS

Complete system consists of (6) "pie-shaped" racks constructed entirely from stainless steel. All systems include choice of durable cardboard, lightweight aluminum, or stainless steel boxes along with choice of 1/2" or 5/8" cell dividers. All systems are available for either liquid phase or vapor phase storage. All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-

PRS - 14 - 9 - 2 - C - 1/2

PRS - Pie Rack System
 PRC - Pie Rack Complete
 PRO - Pie Rack Only

Inside Diameter of Unit into which Racks are fitted

Number of Boxes per Rack

Height of Box

Type of Box: C - Cardboard
 S - Stainless Steel A - Aluminum


Cell Divider Opening 1/2" or 5/8"

COMPLETE SYSTEMS

LIQUID PHASE

VAPOR PHASE

DESCRIPTION	PART NUMBER	DESCRIPTION	PART NUMBER
PRS-14-9-2-C	4101	PRS-14-7-2-C	4105
PRS-14-9-2-S	4102	PRS-14-7-2-S	4106
PRS-14-6-3-C	4103	PRS-14-6-2-C	4107
PRS-14-6-3-S	4104	PRS-14-6-2-S	4108
		PRS-14-4-3-C	4109
		PRS-14-4-3-S	4110





INDIVIDUAL RACKS

COMPLETE

RACKS ONLY

DESCRIPTION	PART NUMBER	DESCRIPTION	PART NUMBER
PRC-14-9-2-C	4201	PRO-14-9-2	4301
PRC-14-9-2-S	4202	PRO-14-7-2	4302
PRC-14-7-2-C	4203	PRO-14-6-2	4303
PRC-14-7-2-S	4204	PRO-14-6-3	4304
PRC-14-6-2-C	4205	PRO-14-4-3	4305
PRC-14-6-2-S	4206		
PRC-14-6-3-C	4207		
PRC-14-6-3-S	4208		
PRC-14-4-3-C	4209		
PRC-14-4-3-S	4210		




ACCESSORIES

INDIVIDUAL BOXES

DESCRIPTION	PART NUMBER
Cardboard - 2"	4B2C
Cardboard - 3"	4B3C
Aluminum - 2"	4B2A
Aluminum - 3"	4B3A
Stainless Steel - 2"	4B2S
Stainless Steel - 3"	4B3S

DISCONTINUED



INDIVIDUAL DIVIDERS

DESCRIPTION	PART NUMBER
1/2" Divider (56 Cell) 12.70 mm opening	4D12
5/8" Divider (42 Cell) 15.87 mm opening	4D58

EXTRAS

DESCRIPTION	PART NUMBER
Platform Riser	4PR
Locking Rods	4LR

NOTES: 1. All vapor phase systems include vapor platforms. 2. Specify divider size when ordering racks.

16" LIQUID NITROGEN FREEZER RACKS

COMPLETE SYSTEMS AND INDIVIDUAL RACKS FOR 16 INCH DIAMETER CONTAINERS

Complete system consists of (6) "pie-shaped" racks constructed entirely from stainless steel. All systems include choice of durable cardboard, lightweight aluminum, or stainless steel boxes along with choice of 1/2" or 5/8" cell dividers. All systems are available for either liquid phase or vapor phase storage. All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-

PRS - 16 - 9 - 2 - C - 1/2

PRS - Pie Rack System
 PRC - Pie Rack Complete
 PRO - Pie Rack Only

Inside Diameter of Unit into which Racks are fitted

Number of Boxes per Rack

Height of Box


Type of Box: C - Cardboard
 S - Stainless Steel A - Aluminum

Cell Divider Opening 1/2" or 5/8"

COMPLETE SYSTEMS

LIQUID PHASE


DESCRIPTION	PART NUMBER	DESCRIPTION	PART NUMBER
PRS-16-12-2-C	6101	PRS-16-9-2-C	6119
PRS-16-12-2-S	6102	PRS-16-9-2-S	6120
PRS-16-9-2-C	6107	PRS-16-6-2-C	6121
PRS-16-9-2-S	6108	PRS-16-6-2-S	6122
PRS-16-8-3-C	6109	PRS-16-6-3-C	6125
PRS-16-8-3-S	6110	PRS-16-6-3-S	6126
PRS-16-6-3-C	6111	PRS-16-4-3-C	6127
PRS-16-6-3-S	6112	PRS-16-4-3-S	6128
PRS-16-6-3.75-C	6113	PRS-16-5-3.75-C	6131
PRS-16-6-3.75-S	6114	PRS-16-5-3.75-S	6132
PRS-16-5-3.75-C	6115	PRS-16-4-3.75-C	6133
PRS-16-5-3.75-S	6116	PRS-16-4-3.75-S	6134



INDIVIDUAL RACKS


COMPLETE

DESCRIPTION	PART NUMBER
PRC-16-12-2-C	6201
PRC-16-12-2-S	6202
PRC-16-9-2-C	6207
PRC-16-9-2-S	6208
PRC-16-6-2-C	6209
PRC-16-6-2-S	6210
PRC-16-8-3-C	6211
PRC-16-8-3-S	6212
PRC-16-6-3-C	6215
PRC-16-6-3-S	6216
PRC-16-4-3-C	6217
PRC-16-4-3-S	6218
PRC-16-6-3.75-C	6219
PRC-16-6-3.75-S	6220
PRC-16-5-3.75-C	6221
PRC-16-5-3.75-S	6222
PRC-16-4-3.75-C	6223
PRC-16-4-3.75-S	6224



RACKS ONLY


DESCRIPTION	PART NUMBER
PRO-16-12-2	6301
PRO-16-9-2	6304
PRO-16-6-2	6305
PRO-16-8-3	6306
PRO-16-6-3	6308
PRO-16-4-3	6309
PRO-16-6-3.75	6310
PRO-16-5-3.75	6311
PRO-16-4-3.75	6312



ACCESSORIES

INDIVIDUAL BOXES

DESCRIPTION	PART NUMBER
Cardboard - 2"	6B2C
Cardboard - 3"	6B3C
Aluminum - 2"	6B2A
Aluminum - 3"	6B3A
Stainless Steel - 2"	6B2S
Stainless Steel - 3"	6B3S



INDIVIDUAL DIVIDERS

DESCRIPTION	PART NUMBER
1/2" Divider (81 Cell) 12.70 mm opening	6D12
5/8" Divider (53 Cell) 15.87 mm opening	6D58

EXTRAS

DESCRIPTION	PART NUMBER
Platform Riser	6PR
Locking Rods	6LR


NOTES: 1. All vapor phase systems include vapor platforms. 2. Specify divider size when ordering racks.

18" LIQUID NITROGEN FREEZER RACKS

COMPLETE SYSTEMS AND INDIVIDUAL RACKS FOR 18 INCH DIAMETER CONTAINERS

Complete system consists of (6) "pie-shaped" racks constructed entirely from stainless steel. All systems include choice of durable cardboard, lightweight aluminum, or stainless steel boxes along with choice of 1/2" or 5/8" cell dividers. All systems are available for either liquid phase or vapor phase storage. All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-




COMPLETE SYSTEMS

LIQUID PHASE

VAPOR PHASE

DESCRIPTION	PART NUMBER	DESCRIPTION	PART NUMBER
PRS-18-12-2-C	8101	PRS-18-10-2-C	8109
PRS-18-12-2-S	8102	PRS-18-10-2-S	8110
PRS-18-11-2-C	8103	PRS-18-9-2-C	8111
PRS-18-11-2-S	8104	PRS-18-9-2-S	8112
PRS-18-8-3-C	8105	PRS-18-7-3-C	8113
PRS-18-8-3-S	8106	PRS-18-7-3-S	8114
PRS-18-6-3.75-C	8107	PRS-18-6-3-C	8115
PRS-18-6-3.75-S	8108	PRS-18-6-3-S	8116
		PRS-18-5-3.75-C	8117
		PRS-18-5-3.75-S	8118




INDIVIDUAL RACKS

COMPLETE

RACKS ONLY


DESCRIPTION	PART NUMBER	DESCRIPTION	PART NUMBER
PRC-18-12-2-C	8201	PRO-18-12-2	8301
PRC-18-12-2-S	8202	PRO-18-11-2	8302
PRC-18-11-2-C	8203	PRO-18-10-2	8303
PRC-18-11-2-S	8204	PRO-18-9-2	8304
PRC-18-10-2-C	8205	PRO-18-8-3	8305
PRC-18-10-2-S	8206	PRO-18-7-3	8306
PRC-18-9-2-C	8207	PRO-18-6-3	8307
PRC-18-9-2-S	8208	PRO-18-6-3.75	8308
PRC-18-8-3-C	8209	PRO-18-5-3.75	8309
PRC-18-8-3-S	8210		
PRC-18-7-3-C	8211		
PRC-18-7-3-S	8212		
PRC-18-6-3-C	8213		
PRC-18-6-3-S	8214		
PRC-18-6-3.75-C	8215		
PRC-18-6-3.75-S	8216		
PRC-18-5-3.75-C	8217		
PRC-18-5-3.75-S	8218		



ACCESSORIES

INDIVIDUAL BOXES

DESCRIPTION	PART NUMBER
Cardboard - 2"	8B2C
Cardboard - 3"	8B3C
Aluminum - 2"	8B2A
Aluminum - 3"	8B3A
Stainless Steel - 2"	8B2S
Stainless Steel - 3"	8B3S



INDIVIDUAL DIVIDERS

DESCRIPTION	PART NUMBER
1/2" Divider (112 Cell) 12.70 mm opening	8D12
5/8" Divider 15.87 mm opening	8D58

EXTRAS

DESCRIPTION	PART NUMBER
Platform Riser	8PR
Locking Rods	8LR


NOTES: 1. All vapor phase systems include vapor platforms. 2. Specify divider size when ordering racks.

19" LIQUID NITROGEN FREEZER RACKS

COMPLETE SYSTEMS AND INDIVIDUAL RACKS FOR 19 INCH DIAMETER CONTAINERS

Complete system consists of (7) "pie-shaped" racks constructed entirely from stainless steel. All systems include choice of durable cardboard, lightweight aluminum, or stainless steel boxes along with choice of 1/2" or 5/8" cell dividers. All systems are available for either liquid phase or vapor phase storage. All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-




COMPLETE SYSTEMS

LIQUID PHASE

DESCRIPTION	PART NUMBER	DESCRIPTION	PART NUMBER
PRS-19-11-2-C	9101	PRS-19-9-2-C	9107
PRS-19-11-2-S	9102	PRS-19-9-2-S	9108
PRS-19-7-3-C	9103	PRS-19-6-3-C	9109
PRS-19-7-3-S	9104	PRS-19-6-3-S	9110

VAPOR PHASE





INDIVIDUAL RACKS

COMPLETE

DESCRIPTION	PART NUMBER
PRC-19-11-2-C	9201
PRC-19-11-2-S	9202
PRC-19-9-2-C	9205
PRC-19-9-2-S	9206
PRC-19-7-3-C	9207
PRC-19-7-3-S	9208
PRC-19-6-3-C	9209
PRC-19-6-3-S	9210

RACKS ONLY

DESCRIPTION	PART NUMBER
PRO-19-11-2	9301
PRO-19-9-2	9303
PRO-19-7-3	9304
PRO-19-6-3	9305




ACCESSORIES

INDIVIDUAL BOXES

DESCRIPTION	PART NUMBER
Cardboard - 2"	9B2C
Cardboard - 3"	9B3C
Aluminum - 2"	9B2A
Aluminum - 3"	9B3A
Stainless Steel - 2"	9B2S
Stainless Steel - 3"	9B3S

DISCONTINUED



INDIVIDUAL DIVIDERS

DESCRIPTION	PART NUMBER
1/2" Divider (86 Cell) 12.70 mm opening (Cardboard box only)	9D12C
5/8" Divider (64 Cell) 15.87 mm opening (Cardboard box only)	9D58C
1/2" Divider (104 Cell) 12.70 mm opening (Stainless steel box only)	9D12S
5/8" Divider (80 Cell) 15.87 mm opening (Stainless steel box only)	9D58S

EXTRAS

DESCRIPTION	PART NUMBER
Platform Riser	9PR
Locking Rods	9LR

NOTES: 1. All vapor phase systems include vapor platforms. 2. Specify divider size when ordering racks.

20" LIQUID NITROGEN FREEZER RACKS

COMPLETE RACK SYSTEMS FOR 20 INCH DIAMETER CONTAINERS

Complete system consists of (7) square racks constructed entirely from stainless steel. All systems include choice of durable cardboard, lightweight aluminum, stainless steel or plastic boxes (*soon available*). Each system also includes choice of cell dividers; 100, 81, 64, 49, 36, 25. All systems are available for either liquid phase or vapor phase storage.

All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-

SRS - 20 - 12 - 2 - C - 81

SRS – Square Rack System

Inside Diameter of Unit into which Racks are fitted

Number of Boxes per Rack

Height of Box

Type of Box:
C - Cardboard S - Stainless Steel
A - Aluminum P - Plastic (*soon available*)


Number of Cell Openings:
100, 81, 64, 49, 36, 25, etc.

LIQUID PHASE SYSTEMS

VAPOR PHASE SYSTEMS

DESCRIPTION	PART NUMBER	SYSTEM HEIGHT
SRS-20-13-2-C	2001	27½"
SRS-20-13-2-S	2002	27½"
SRS-20-12-2-C	2003	25½"
SRS-20-12-2-S	2004	25½"
SRS-20-11-2-C	2005	23¾"
SRS-20-11-2-S	2006	23¾"
SRS-20-9-3-C	2007	28⅛"
SRS-20-9-3-S	2008	28⅛"
SRS-20-8-3-C	2009	25"
SRS-20-8-3-S	2010	25"
SRS-20-7-3.75-C	2011	26¾"
SRS-20-7-3.75-S	2012	26¾"
SRS-20-6-3.75-C	2013	23"
SRS-20-6-3.75-S	2014	23"

DESCRIPTION	PART NUMBER	SYSTEM HEIGHT
SRS-20-10-2-C	2015	25½"
SRS-20-10-2-S	2016	25½"
SRS-20-9-2-C	2017	25¼"
SRS-20-9-2-S	2018	25¼"
SRS-20-7-3-C	2019	26"
SRS-20-7-3-S	2020	26"
SRS-20-6-3-C	2021	25"
SRS-20-6-3-S	2022	25"
SRS-20-5-3.75-C	2023	25¼"
SRS-20-5-3.75-S	2024	25¼"



NOTES: 1. All vapor phase systems include vapor platforms. 2. Specify divider size when ordering racks.

31" LIQUID NITROGEN FREEZER RACKS

COMPLETE RACK SYSTEMS FOR 31 INCH DIAMETER CONTAINERS

Complete system consists of (17) square racks constructed entirely from stainless steel. All systems include choice of durable cardboard, lightweight aluminum, stainless steel or plastic boxes (*soon available*). Each system also includes choice of cell dividers; 100, 81, 64, 49, 36, 25. All systems are available for either liquid phase or vapor phase storage.

All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-

SRS - 31 - 12 - 2 - C - 81

SRS – Square Rack System

Inside Diameter of Unit into which Racks are fitted

Number of Boxes per Rack

Height of Box

Type of Box:
C - Cardboard S - Stainless Steel
A - Aluminum P - Plastic (soon available)


Number of Cell Openings:
100, 81, 64, 49, 36, 25, etc.

LIQUID PHASE SYSTEMS

DESCRIPTION	PART NUMBER	SYSTEM HEIGHT
SRS-31-13-2-C	3101	27½"
SRS-31-13-2-S	3102	27½"
SRS-31-12-2-C	3103	25½"
SRS-31-12-2-S	3104	25½"
SRS-31-11-2-C	3105	23 ⁷ / ₁₆ "
SRS-31-11-2-S	3106	23 ⁷ / ₁₆ "
SRS-31-9-3-C	3107	28 ¹ / ₈ "
SRS-31-9-3-S	3108	28 ¹ / ₈ "
SRS-31-8-3-C	3109	25"
SRS-31-8-3-S	3110	25"
SRS-31-7-3.75-C	3111	26 ³ / ₄ "
SRS-31-7-3.75-S	3112	26 ³ / ₄ "
SRS-31-6-3.75-C	3113	23"
SRS-31-6-3.75-S	3114	23"

VAPOR PHASE SYSTEMS

DESCRIPTION	PART NUMBER	SYSTEM HEIGHT
SRS-31-10-2-C	3115	25½"
SRS-31-10-2-S	3116	25½"
SRS-31-9-2-C	3117	25¼"
SRS-31-9-2-S	3118	25¼"
SRS-31-7-3-C	3119	26"
SRS-31-7-3-S	3120	26"
SRS-31-6-3-C	3121	25"
SRS-31-6-3-S	3122	25"
SRS-31-5-3.75-C	3123	25¼"
SRS-31-5-3.75-S	3124	25¼"



NOTES: 1. All vapor phase systems include vapor platforms.

2. Specify divider size when ordering racks.

39" LIQUID NITROGEN FREEZER RACKS

COMPLETE RACK SYSTEMS FOR 39 INCH DIAMETER CONTAINERS

Complete system consists of (28) square racks constructed entirely from stainless steel. All systems include choice of durable cardboard, lightweight aluminum, stainless steel or plastic boxes (*soon available*). Each system also includes choice of cell dividers; 100, 81, 64, 49, 36, 25. All systems are available for either liquid phase or vapor phase storage.

All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-

SRS - 39 - 12 - 2 - C - 81

SRS - Square Rack System

Inside Diameter of Unit into which Racks are fitted

Number of Boxes per Rack

Height of Box

Type of Box:

C - Cardboard S - Stainless Steel

A - Aluminum P - Plastic (soon available)


Number of Cell Openings: 100, 81, 64, 49, 36, 25, etc.

LIQUID PHASE SYSTEMS

DESCRIPTION	PART NUMBER	SYSTEM HEIGHT
SRS-39-13-2-C	3301	27½"
SRS-39-13-2-S	3302	27½"
SRS-39-12-2-C	3303	25½"
SRS-39-12-2-S	3304	25½"
SRS-39-11-2-C	3305	23¾"
SRS-39-11-2-S	3306	23¾"
SRS-39-9-3-C	3307	28⅛"
SRS-39-9-3-S	3308	28⅛"
SRS-39-8-3-C	3309	25"
SRS-39-8-3-S	3310	25"
SRS-39-7-3.75-C	3311	26¾"
SRS-39-7-3.75-S	3312	26¾"
SRS-39-6-3.75-C	3313	23"
SRS-39-6-3.75-S	3314	23"

VAPOR PHASE SYSTEMS

DESCRIPTION	PART NUMBER	SYSTEM HEIGHT
SRS-39-10-2-C	3315	25½"
SRS-39-10-2-S	3316	25½"
SRS-39-9-2-C	3317	25¼"
SRS-39-9-2-S	3318	25¼"
SRS-39-7-3-C	3319	26"
SRS-39-7-3-S	3320	26"
SRS-39-6-3-C	3321	25"
SRS-39-6-3-S	3322	25"
SRS-39-5-3.75-C	3323	25¼"
SRS-39-5-3.75-S	3324	25¼"




NOTES: 1. All vapor phase systems include vapor platforms.

2. Specify divider size when ordering racks.

LIQUID NITROGEN FREEZER RACKS

INDIVIDUAL SQUARE RACKS



SRO-8-3

SRO-12-2

SRO-5-3.75

Individual Square Racks are constructed entirely from stainless steel. All racks are available complete with boxes and cell dividers or purchased separately. All square racks have a standard 5 1/2" x 5 1/2" (front to back), (side to side) dimension. All racks include locking rods. Custom racks available upon request.

-EXAMPLE OF DESCRIPTION-

SRC - 12 - 2 - C - 81

SRC - Square Rack Complete
 SRO - Square Rack Only
 MRO - Mini Rack Only

Number of Boxes per Rack: 12
 Height of Box: 2
 Type of Box: C - Cardboard
 S - Stainless Steel A - Aluminum

Number of Cell Openings: 100, 81, 64, 49, 25, 16

STANDARD RACKS COMPLETE

DESCRIPTION	PART NUMBER	CAPACITY	RACK HEIGHT	
			STANDARD	OPTIONAL
SRC-13-2-C	1101	(13) Boxes	27 1/2"	28 1/2"
SRC-13-2-S	1102	(13) Boxes	27 1/2"	28 1/2"
SRC-12-2-C	1103	(12) Boxes	25 1/2"	26 9/16"
SRC-12-2-S	1104	(12) Boxes	25 1/2"	26 9/16"
SRC-11-2-C	1105	(11) Boxes	23 7/16"	24 1/2"
SRC-11-2-S	1106	(11) Boxes	23 7/16"	24 1/2"
SRC-10-2-C	1107	(10) Boxes	21 3/8"	22 1/4"
SRC-10-2-S	1108	(10) Boxes	21 3/8"	22 1/4"
SRC-9-2-C	1109	(9) Boxes	19 1/4"	20"
SRC-9-2-S	1110	(9) Boxes	19 1/4"	20"
SRC-8-2-C	1111	(8) Boxes	17 1/8"	17 3/4"
SRC-8-2-S	1112	(8) Boxes	17 1/8"	17 3/4"
SRC-7-2-C	1113	(7) Boxes	15"	15 5/8"
SRC-7-2-S	1114	(7) Boxes	15"	15 5/8"
SRC-6-2-C	1115	(6) Boxes	12 15/16"	13 1/16"
SRC-6-2-S	1116	(6) Boxes	12 15/16"	13 1/16"
SRC-5-2-C	1117	(5) Boxes	10 3/4"	11 5/16"
SRC-5-2-S	1118	(5) Boxes	10 3/4"	11 5/16"
SRC-4-2-C	1119	(4) Boxes	8 1/16"	9 1/8"
SRC-4-2-S	1120	(4) Boxes	8 1/16"	9 1/8"
SRC-3-2-C	1121	(3) Boxes	6 5/8"	7"
SRC-3-2-S	1122	(3) Boxes	6 5/8"	7"
SRC-9-3-C	1123	(9) Boxes	28 1/8"	28 7/8"
SRC-9-3-S	1124	(9) Boxes	28 1/8"	28 7/8"
SRC-8-3-C	1125	(8) Boxes	25"	25 7/8"
SRC-8-3-S	1126	(8) Boxes	25"	25 7/8"
SRC-7-3-C	1127	(7) Boxes	21 15/16"	22 3/4"
SRC-7-3-S	1128	(7) Boxes	21 15/16"	22 3/4"
SRC-6-3-C	1129	(6) Boxes	18 7/8"	19 1/2"
SRC-6-3-S	1130	(6) Boxes	18 7/8"	19 1/2"
SRC-5-3-C	1131	(5) Boxes	15 3/4"	16 3/8"
SRC-5-3-S	1132	(5) Boxes	15 3/4"	16 3/8"
SRC-4-3-C	1133	(4) Boxes	12 5/8"	13 1/4"
SRC-4-3-S	1134	(4) Boxes	12 5/8"	13 1/4"
SRC-7-3.75-C	1135	(7) Boxes	26 3/4"	26 3/4"
SRC-7-3.75-S	1136	(7) Boxes	26 3/4"	26 3/4"
SRC-6-3.75-C	1137	(6) Boxes	23"	23"
SRC-6-3.75-S	1138	(6) Boxes	23"	23"
SRC-5-3.75-C	1139	(5) Boxes	19 1/8"	19 1/8"
SRC-5-3.75-S	1140	(5) Boxes	19 1/8"	19 1/8"
SRC-4-3.75-C	1141	(4) Boxes	15 3/8"	15 3/8"
SRC-4-3.75-S	1142	(4) Boxes	15 3/8"	15 3/8"

STANDARD RACKS ONLY

DESCRIPTION	PART NUMBER	CAPACITY	RACK HEIGHT	
			STANDARD	OPTIONAL
SRO-13-2	1201	(13) Boxes	27 1/2"	28 1/2"
SRO-12-2	1202	(12) Boxes	25 1/2"	26 9/16"
SRO-11-2	1203	(11) Boxes	23 7/16"	24 1/2"
SRO-10-2	1204	(10) Boxes	21 3/8"	22 1/4"
SRO-9-2	1205	(9) Boxes	19 1/4"	20"
SRO-8-2	1206	(8) Boxes	17 1/8"	17 3/4"
SRO-7-2	1207	(7) Boxes	15"	15 5/8"
SRO-6-2	1208	(6) Boxes	12 15/16"	13 1/16"
SRO-5-2	1209	(5) Boxes	10 3/4"	11 5/16"
SRO-4-2	1210	(4) Boxes	8 1/16"	9 1/8"
SRO-3-2	1211	(3) Boxes	6 5/8"	7"
SRO-9-3	1212	(9) Boxes	28 1/8"	28 7/8"
SRO-8-3	1213	(8) Boxes	25"	25 7/8"
SRO-7-3	1214	(7) Boxes	21 15/16"	22 3/4"
SRO-6-3	1215	(6) Boxes	18 7/8"	19 1/2"
SRO-5-3	1216	(5) Boxes	15 3/4"	16 3/8"
SRO-4-3	1217	(4) Boxes	12 5/8"	13 1/4"
SRO-7-3.75	1218	(7) Boxes	26 3/4"	26 3/4"
SRO-6-3.75	1219	(6) Boxes	23"	23"
SRO-5-3.75	1220	(5) Boxes	19 1/8"	19 1/8"
SRO-4-3.75	1221	(4) Boxes	15 3/8"	15 3/8"

MINI RACKS ONLY

DESCRIPTION	PART NUMBER	CAPACITY	RACK HEIGHT	
			STANDARD	OPTIONAL
MRO-13-2	M101	(13) Boxes	27 1/2"	28 1/2"
MRO-12-2	M102	(12) Boxes	25 1/2"	26 9/16"
MRO-11-2	M103	(11) Boxes	23 7/16"	24 1/2"
MRO-10-2	M104	(10) Boxes	21 3/8"	22 1/4"
MRO-9-2	M105	(9) Boxes	19 1/4"	20"
MRO-8-2	M106	(8) Boxes	17 1/8"	17 3/4"
MRO-7-2	M107	(7) Boxes	15"	15 5/8"
MRO-6-2	M108	(6) Boxes	12 15/16"	13 1/16"
MRO-5-2	M109	(5) Boxes	10 3/4"	11 5/16"
MRO-4-2	M110	(4) Boxes	8 1/16"	9 1/8"
MRO-3-2	M111	(3) Boxes	6 5/8"	7"
MRO-9-3	M112	(9) Boxes	28 1/8"	28 7/8"
MRO-8-3	M113	(8) Boxes	25"	25 7/8"
MRO-7-3	M114	(7) Boxes	21 15/16"	22 3/4"
MRO-6-3	M115	(6) Boxes	18 7/8"	19 1/2"
MRO-5-3	M116	(5) Boxes	15 3/4"	16 3/8"
MRO-4-3	M117	(4) Boxes	12 5/8"	13 1/4"
MRO-7-3.75	M118	(7) Boxes	26 3/4"	26 3/4"
MRO-6-3.75	M119	(6) Boxes	23"	23"
MRO-5-3.75	M120	(5) Boxes	19 1/8"	19 1/8"
MRO-4-3.75	M121	(4) Boxes	15 3/8"	15 3/8"

NOTES: 1. Specify divider size when ordering racks. 2. Rack height (optional), accepts plastic boxes also.

STANDARD BOXES AND CELL DIVIDERS


For All Upright Freezer And Chest Freezer Racks.

Constructed of recycled paper products, aluminum or stainless steel. Custom sizes available upon request.

STANDARD BOXES


CARDBOARD

<u>DESCRIPTION</u>	<u>PART NUMBER</u>	<u>DIMENSIONS (OUTSIDE)</u>
Regular - 2"	B2C	5¼" x 5¼" x 1⅞"
Regular - 3"	B3C	5¼" x 5¼" x 2⅞"




STAINLESS STEEL

<u>DESCRIPTION</u>	<u>PART NUMBER</u>	<u>DIMENSIONS (OUTSIDE)</u>
Regular - 2"	B2S	5" x 5" x 1⅞"
Regular - 3"	B3S	5" x 5" x 2⅞"



ALUMINUM


<u>DESCRIPTION</u>	<u>PART NUMBER</u>	<u>DIMENSIONS (OUTSIDE)</u>
Regular - 2"	B2A	5" x 5" x 1⅞"
Regular - 3"	B3A	5" x 5" x 2⅞"



STANDARD CELL DIVIDERS


100 CELL

<u>PART NUMBER</u>	<u>CELL OPENING</u>	<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u> <u>IN</u>	
D100	11.65 0.45	4⅞" x 4⅞"




81 CELL

<u>PART NUMBER</u>	<u>CELL OPENING</u>	<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u> <u>IN</u>	
D81	12.70 0.50	4⅞" x 4⅞"




64 CELL

<u>PART NUMBER</u>	<u>CELL OPENING</u>	<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u> <u>IN</u>	
D64	14.60 0.57	4⅞" x 4⅞"




49 CELL

<u>PART NUMBER</u>	<u>CELL OPENING</u>	<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u> <u>IN</u>	
D49	16.76 0.66	4⅞" x 4⅞"




25 CELL

<u>PART NUMBER</u>	<u>CELL OPENING</u>	<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u> <u>IN</u>	
D25	24.51 0.96	4⅞" x 4⅞"



16 CELL

<u>PART NUMBER</u>	<u>CELL OPENING</u>	<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u> <u>IN</u>	
D16	30.23 1.19	4⅞" x 4⅞"



NOTE: Metal boxes available with lift off or hinged tops. Consult factory for prices.

OPTIONAL BOXES AND CELL DIVIDERS


For All Upright Freezer And Chest Freezer Racks.

Constructed of recycled paper products, aluminum or stainless steel. Custom sizes available upon request.

OPTIONAL BOXES


CARDBOARD

<u>DESCRIPTION</u>	<u>PART NUMBER</u>	<u>DIMENSIONS (OUTSIDE)</u>
Large - 2"	B2CL	5 $\frac{3}{4}$ " x 5 $\frac{3}{4}$ " x 1 $\frac{7}{8}$ "
Large - 3"	B3CL	5 $\frac{3}{4}$ " x 5 $\frac{3}{4}$ " x 2 $\frac{7}{8}$ "
Special - 2"	B2CS	5 $\frac{5}{8}$ " x 5 $\frac{5}{8}$ " x 1 $\frac{7}{8}$ "
Special - 3"	B3CS	5 $\frac{5}{8}$ " x 5 $\frac{5}{8}$ " x 2 $\frac{7}{8}$ "
Mini - 2"	B2CM	3" x 3" 1 $\frac{7}{8}$ "
Mini - 3"	B3CM	3" x 3" x 2 $\frac{7}{8}$ "




STAINLESS STEEL

<u>DESCRIPTION</u>	<u>PART NUMBER</u>	<u>DIMENSIONS (OUTSIDE)</u>
Large - 2"	B2SL	5 $\frac{5}{8}$ " x 5 $\frac{5}{8}$ " x 1 $\frac{7}{8}$ "
Large - 3"	B3SL	5 $\frac{5}{8}$ " x 5 $\frac{5}{8}$ " x 2 $\frac{7}{8}$ "
Special - 2"	B2SS	5 $\frac{1}{4}$ " x 5 $\frac{1}{4}$ " x 1 $\frac{7}{8}$ "
Special - 3"	B3SS	5 $\frac{1}{4}$ " x 5 $\frac{1}{4}$ " x 2 $\frac{7}{8}$ "
Mini - 2"	B2SM	3" x 3" x 1 $\frac{7}{8}$ "
Mini - 3"	B3SM	3" x 3" x 2 $\frac{7}{8}$ "



ALUMINUM


<u>DESCRIPTION</u>	<u>PART NUMBER</u>	<u>DIMENSIONS (OUTSIDE)</u>
Large - 2"	B2AL	5 $\frac{5}{8}$ " x 5 $\frac{5}{8}$ " x 1 $\frac{7}{8}$ "
Large - 3"	B3AL	5 $\frac{5}{8}$ " x 5 $\frac{5}{8}$ " x 2 $\frac{7}{8}$ "
Special - 2"	B2AS	5 $\frac{1}{4}$ " x 5 $\frac{1}{4}$ " x 1 $\frac{7}{8}$ "
Special - 3"	B3AS	5 $\frac{1}{4}$ " x 5 $\frac{1}{4}$ " x 2 $\frac{7}{8}$ "
Mini - 2"	B2AM	3" x 3" x 1 $\frac{7}{8}$ "
Mini - 3"	B3AM	3" x 3" x 2 $\frac{7}{8}$ "



OPTIONAL CELL DIVIDERS


100 CELL LARGE

<u>PART NUMBER</u>	<u>CELL OPENING</u>		<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u>	<u>IN</u>	
D100L	12.70	0.50	5 $\frac{1}{2}$ " x 5 $\frac{1}{2}$ "




100 CELL SPECIAL

<u>PART NUMBER</u>	<u>CELL OPENING</u>		<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u>	<u>IN</u>	
D100S	12.70	0.50	5 $\frac{1}{8}$ " x 5 $\frac{1}{8}$ "



25 CELL MINI

<u>PART NUMBER</u>	<u>CELL OPENING</u>		<u>DIMENSIONS (OUTSIDE)</u>
	<u>MM</u>	<u>IN</u>	
D25M	12.70	0.50	2 $\frac{7}{8}$ " x 2 $\frac{7}{8}$ "



NOTE: Metal boxes available with lift off or hinged tops. Consult factory for prices.



Viewing and Printing Tips for Acrobat PDF Documents

Note: the interface and specific techniques may vary depending on the version of Acrobat Reader that you are using. Adobe Acrobat Reader is a free program. You may download the latest version from Adobe.com. To see complete instructions for Acrobat, use the HELP button at the top of Acrobat.

BETTER VIEWING

The clarity of a PDF document on your screen may vary depending on the size you choose to view the document at.

You can read dimensions better and generally view a document better by zooming in. To zoom in, click the Zoom In button on the top toolbar (the + sign). Or, choose Zoom In from the View menu. You can decrease the size of the document with the Zoom Out tool (the - sign)

To zoom in on a specific area, click the Zoom button on the toolbar, then click in the area on the document which you want to enlarge. Zooming out works the same way.

You can also view a document at a specific zoom level. Select Zoom To from the View menu. The zoom dialog box appears. Type the zoom percentage you desire in the text box and click OK. The document will be displayed at your specified magnification.

PRINTING

Regardless of the magnification and clarity of the document you are viewing, printing the document will generally result in a very clear image of text and product dimensions. These Acrobat documents are optimized to print clearly on your local computer printer.

To print the document to your local computer printer, either choose Print from the File menu, click the Print button on the File toolbar, or press Ctrl P (Command P on the Mac). The Print dialog box appears.

You may select the number of copies you want to print in the Print dialog box.

If you don't want to print the entire document, enter the page range (the starting and ending page numbers of the pages to print) in the appropriate text boxes of the Print dialog box.

Click OK (Windows) or Print (Mac)

For additional printing instructions, use the Help button on the top toolbar of Acrobat Reader.

Note: if you need to print this catalog in part or in its entirety through a commercial printer, contact Custom BioGenic Systems for an Adobe Acrobat PDF version suitable for high-end commercial printing.