

PORTALL™ VERSION 2

Software for particle counters



Schedule and automate handheld and portable devices

Features

- Supports 21 CFR Part 11 compliance
- Intuitive Windows®-based software dramatically reduces operator training time
- Track particle count data, identify trends and eliminate particle contamination sources

Working with your handheld and portable device, PORTALL VERSION 2 simply and easily organizes archives and retrieves particle count data.

Gown up less

PortAll's scheduling feature automates particle counting and data collection. Get inside data conveniently from outside the cleanroom.

Report with ease

A Windows-based software interface and report wizard makes calculation and presentation of ISO14644 and FS209E easy. Name individual particle counters to provide reporting from multiple counters for the same location.

Simplify validation

Verify adherence to standards and regulations, including ISO and 21 CFR Part 11. The software also supports conducting rapid pass/fail reporting and documentation.

Expand your particle counter's capabilities

PORTALL VERSION 2 expands the functionality of PortAll 1.0 to include:

- Communication with Met One HHPG-2 and HHPG-6 particle counters
- Password authentication and management
- Auto logoff
- Lockout after failed log in attempts
- Secure data
- BS5295 calculation
- Pre-configuration of particle counters
- User-defined counter identifier

Performance Specifications

Supported Particle Counters	Handhelds	Portables	Portables for Water
	HHPC-2	237	WGS267
	HHPC-6	2100/2200	
	227	2400/2408	
		3313/3315	
		1104 CNC	
Operating System Software	Microsoft® Windows® 2000 Professional (SP1 or greater) Microsoft Windows XP® Professional Microsoft Internet Explorer® 5.x or greater		
Hardware Minimum Recommended	CPU type Intel® Pentium® III, 500 MHz Intel Pentium IV, 1 GHz Free Disk Space 500 MB 4 GB Memory / RAM 128 MB 256 MB Video Display 800x600* 1024x768 Serial/COM ports 1 RS-232 2 RS-232 or 2 USB		

* Non-English Windows systems require a minimum video display of 1024x768

Global Headquarters

6, route de Compois - CP 212
1222 Vérenaz - Geneva - Switzerland
Tel ++ 41 (0)22 594 64 00
Fax ++ 41 (0)22 594 64 99

Americas Headquarters

481 California Avenue
Grants Pass - Oregon 97526 - USA
Tel 1 800 866 7889 / 1 541 472 6500
Fax 1 541 472 6170



LANGE 

UNITED FOR WATER QUALITY

