

# Pharmaceutical Freezer-MED-70

With temperature indicator • stepless temperature selection



Capacity	70 Litres
Temperature setting	approx. -15° to - 22°C
Voltage	220-240 V, 50/60 Hz
Power draft	100 Watt
Average consumption	1,98 kWh/24 h
Overall dimensions incl. distance from wall	w x d x h = 47 x 51 x 74 cm
Interior dimensions	w x d x h = 35 x 41 x 55 cm (usable depth on top 9 cm, below less 20 cm)
Overall dimensions with door open at 90°	w x d = 47 x 96 cm
Shelf dimensions	20 x 35 cm resp. 31 x 35 cm
Drawer inner dimensions	w x d x h = 29 x 26 x 5 cm
Max. loading capacity per drawer/shelf	9,5 kg/25 kg
Weight	net 34 kg, gross 42 kg

**Outer panels** of galvanised sheet steel (protection against sub surface rusting, white, powder coated).

**Interior** of impact resistant white plastic with moulded in shelf guides 3 shelves, 1 drawer.

**Insulating**, 55 mm, of high-grade pressure foamed and pro-environmental material, extra thick for energy saving effect.

**Door** with plastic magnetic sealing, lockable. Door does not project at side when opened.

**Door hinge** changeable (also subsequently), standard version with hinges on right (as illustrated), deliverable also with door hinge on left at no extra charge.

**Defrosting** manual. Melt water for manual emptying.

**Temperature gear** by interior thermostat. The temperature is kept automatically, independent of changes in ambient temperatures.

**Compressor** static ventilated, hermetically enclosed, energy saving, low noise, easy to service, for 220-240 V AC, 50/60 Hz mains. Other voltages on request. Length of mains cable supplied approx. 1,8 m.

## Equipment at extra charge:

- Pen-recording thermometer for placing loosely inside appliance with 2 waxed paper strips
- Spare waxed paper strips in packs of 10 rolls