

GLASSWARE WASHERS GW1050



GW1050 TECHNICAL FEATURES

STANDARD PROGRAMS INSTALLED	7 (PLUS 1 SERVICE PROGRAM)
MICROPROCESSOR CAPACITY	8 KBYTE
ALPHANUMERICAL DISPLAY	YES
PHASE PARAMETERS MEMORISED	PRESET
WATER TEMPERATURE	MAX 95°C
DOSING PUMPS FOR NEUTRALIZING AGENT (0-50 ml)	YES
DOSING PUMPS FOR ALKALINE DETERGENT (0-50 ml)	OPTIONAL
SAFETY DEVICES	DOOR OPENING, PUMP AND HEATER FAILURE, WATER LEVEL
ALARM INDICATIONS	9
WATER SUPPLY (PRESSURE 1-3.5 BAR)	
COLD MAINS WATER HARDNESS	MAX 42°F
DEMINERALIZED WATER CONDUCTIVITY	<20 mS/cm
BOOSTER PUMP FOR DEMINERALIZED WATER	OPTIONAL
BUILT-IN WATER SOFTENER	YES
RECIRCULATION PUMP	400 L/MIN
WATER HEATER	
ELECTRIC	6 kW MAX
DIMENSIONS LxDxH (mm)	
EXTERNAL (WITH TOP FOR UNDERCOUNTER)	600x600x850 (830)
INTERNAL	520x515x545
NET WEIGHT (kg)	65
STEEL	
WASH TANK	AISI 316L
EXTERNAL PANELS	AISI 304
POWER SUPPLY	
VOLTAGE / MAX POWER	GW1050-1 1/N/PE 230 V – 50 Hz 3.3 kW GW1050 3/N/PE 400 V – 50 Hz 7.0 kW
NOISE LEVEL	50dBA
CONFORMITY TO STANDARDS	72/23/EEC, 93/68/EEC, 89/336/EEC EN61010-1, EN61010-2-045, EN61326:1997+A1:1998

