

## ANIMAL CAGE AND FEEDING BOTTLE WASHER AC 1200

*Steelco AC 1200 is an efficient system able to wash and disinfect on two levels all types of feeding bottles and animal cages, granting an easy maintenance of bottles and cages cleanliness and effective cleaning results with the thermal/chemical destruction of micro-organisms by rinsing at 85 °C through an independent rinsing circuit.*

### THE STRENGTHS OF STEELCO AC 1200:

- × chamber entirely made in stainless steel AISI 316 L, external body AISI 304.
- × construction system with total removal of angle and breeding ground dirty
- × self cleaning tank with rounded edges
- × water filtering system for an excellent washing result
- × two rotating washing and rinsing arms positioned both at the top and at the bottom of the washing chamber plus other two in the middle level basket. The AC 1200 can be used for the treatment of feeding bottles which can be positioned on dedicated supports for a separate rinsing with injectors directly into each bottle (3 supports each of them for 18 bottles for each level can be accepted)
- × manual single door and safety device
- × door made in stainless steel AISI 316
- × disinfection by rinsing at 85 °C
- × rinsing tank capacity 35 lt
- × tank temperature check
- × check pressure device of the washing circuit
- × equipped with 2 dosing pumps of chemical products, under a complete micro-processor control, with the possibility to modify the required products quantity in each cycle
- × flow meter for checking and validation of chemical product quantity
- × check of chemical product quantity, with minimum level alarm
- × RS 232 port provided with printer connection for monitoring and validating washing phases.
- × 40 washing programmes available
- × adjustable and personalized cycles (complete cycle in 5 minutes)
- × exhaust air connection diam. 120 mm

### OPTIONALS:

- × 3rd chemical dosing pump
- × Printer
- × Steam heating

REV.01

## TECHNICAL DATA:

External dimensions ..... mm 1154x1000x1963h  
External dimensions with open doors ..... mm 1154x1469x2295h  
Washing chamber dimensions ..... mm 1550x900x815h  
Electrical connection ..... 400Vac 3N 50Hz  
(on request special electrical connection available)

Tank heating elements kW 18;  
Steam heated - Optional;  
Washing circuit with 1 pump of 2 kW;  
40 lt. rinsing tank 18 kW, provided with pressure pump;  
Total power 20,5 kW;  
Water connection: cold, warm soft water and DI water.

REV.01