


Accessories

Fixed height loading carriage and loading car with 2 shelves

The loading car with 2 shelves are made of 316 stainless steel. The shelf supports are adjustable for load height adjustment (see quote P1). Shelves can hold 100 Kg (220 lbs) evenly distributed on the shelf.

The fixed height loading carriage, made of 304 stainless steel, is provided with a safety lock that avoids the car to accidentally fall off; the same system is used to lock the loading carriage to the sterilizer.



Model	A	B	C1	C2	D	E	Shelves useful loading dimensions			H2	P1
							H1	K	L		
VS 161626 L	531 21	953 37 1/2	1300 51 1/6	1011 39 4/5	871 34 2/7	335 13 1/5	347 13 2/3	341 13 3/7	577 22 5/7	94 3 5/7	40 1 4/7
VS 202038 L	627 24 2/3	904 35 3/5	1354 53 1/3	960 37 4/5	1176 46 2/7	335 13 1/5	448 17 2/3	442 17 2/5	882 34 5/7	94 3 5/7	40 1 4/7
VS 262639 L	627 24 2/3	904 35 3/5	1504 59 1/5	960 37 4/5	1176 46 2/7	488 19 1/5	600 23 5/8	594 23 2/5	922 36 2/7	124 4 7/8	40 1 4/7
VS 262649 L	928 36 1/2	904 35 3/5	1506 59 2/7	960 37 4/5	1478 58 1/5	495 19 1/2	600 23 5/8	594 23 2/5	1177 46 1/3	124 4 7/8	40 1 4/7
VS 262661 L	1233 48 1/2	904 35 3/5	1506 59 2/7	960 37 4/5	1783 70 1/5	495 19 1/2	600 23 5/8	594 23 2/5	1482 58 1/3	124 4 7/8	40 1 4/7


Additional fixed height loading carriage for unload side

Upon request is possible to order an additional fixed height loading carriage for handling of loads between sterilizer and sterile storage.

Perforated shelves in lieu of car and carriage system

The fixed rack with perforated shelf supports is mounted on each side of the sterilizer chamber walls. The sliding shelf supports are adjustable for load height adjustment (see quote P2). Rack and shelves are made of 316 stainless steel. Shelves can hold 100 Kg (220 lbs) evenly distributed on the shelf.

A-A



Model	Shelves useful loading dimensions						P2
	F	G1	G2	H3	H4	H5	
VS 161626 L	533 21	295 11 3/5	261 10 2/7	158 6 2/9	201 8	369 14 1/2	55 2 1/6
VS 202038 L	838 33	398 15 2/3	366 14 2/5	209 8 1/5	252 10	471 18 1/2	70 2 3/4
VS 262639 L	863 34	548 21 4/7	519 20 3/7	285 11 2/9	328 13	623 24 1/2	70 2 3/4
VS 262649 L	1117 44	548 21 4/7	519 20 3/7	285 11 2/9	328 13	623 24 1/2	70 2 3/4
VS 262661 L	1422 56	548 21 4/7	519 20 3/7	285 11 2/9	328 13	623 24 1/2	70 2 3/4